Candidate Name:		Class Adm No
	millennia institute	Shift Laboratory

2021 Preliminary ExaminationsPre-University 3

H2 CHEMISTRY

9729/04

Paper 4 Practical

1 Sept 2021

2 hours 30 min

Candidates answer on the Question paper.

READ THESE INSTRUCTIONS FIRST

Do not turn over this question paper until you are told to do so

Write your name, class and admission number on all the work you hand in.

Write in dark blue or black pen.

You may use an HB pencil for any diagrams or graphs.

Do not use staples, paper clips, glue or correction fluid.

Answer all questions in the spaces provided on the Question Paper.

The use of an approved scientific calculator is expected, where appropriate.

You may lose marks if you do not show your working or if you do not use appropriate units.

Qualitative Analysis Notes are printed at the back of the Question Paper.

At the end of the examination, fasten all your work securely together.

The number of marks is given in brackets [] at the end of each question or part question.

Question	1	2	3	4	Total
Marks	13	20	10	12	55

1 Determine the formula of the ion IO_x by titration

In this experiment you will determine the formula of the ion, IOx-.

lodine, I₂, is formed when IO_x-ions reacts with an excess of iodide ions.

The equation for this reaction is:

$$IO_x^- + y I^- + z H^+ \rightarrow (\frac{y+1}{2}) I_2 + \frac{z}{2} H_2O$$

where x, y and z are all integers.

The amount of iodine produced will then be determined by titration with thiosulfate ions, S₂O₃²⁻.

$$I_2 + 2S_2O_3^{2-} \rightarrow 2I^- + S_4O_6^{2-}$$

FA 1 is a solution containing $0.0150 \text{ mol dm}^{-3} \text{ IO}_x^-$ ions.

FA 2 is dilute sulfuric acid, H₂SO₄.

FA 3 is 0.500 mol dm⁻³ potassium iodide, KI.

FA 4 is 0.100 mol dm⁻³ sodium thiosulfate, Na₂S₂O₃.

starch indicator

(a) Method

- Pipette 25.0 cm³ of FA 1 into a conical flask.
- Using a measuring cylinder, add 20.0 cm³ of FA 2 to the conical flask.
- Using a measuring cylinder, add 10.0 cm³ of FA 3 to the conical flask. The solution will turn brown as iodine is produced.
- Fill the burette with FA 4.
- Add FA 4 from the burette until the solution in the conical flask turns yellow.
- Add 1 cm³ of starch indicator to the conical flask. The solution will turn blue-black.
- Continue to add more FA 4 dropwise from the burette until the blue-black colour just disappears.

This is the end-point of the titration.

- Carry out as many accurate titrations as you think necessary to obtain consistent results.
- Record the burette readings and the volume of FA 4 added in each accurate titration in the space below.

Results

M1

M2

М3

Μ4

M5

(b)	From your titration results, obtain a value for the volume of FA 4 to be used in your	
	calculations. Show clearly how you obtained this value.	
	25.0 cm ³ of FA 1 required cm ³ of FA 4 .	M6
(c)	Calculations	
(0)		
	(i) Use your answer to (b) to calculate the amount of iodine that was formed.	
	·	
	Amount of $I_2 = \dots$	M7
	, and an ing the second of the	
	(ii) Calculate the amount of IO _x ⁻ ions in 25.0 cm ³ of FA 1.	
	(ii) Calculate the amount of TO_X (one ii) 20.0 cm (STA).	
		ĺ
		M8
	amount of IO _x ⁻ ions=	
		ļ

	(iii)	Using your answers in (c)(i) and (c)(ii), calculate the value of y.	
		$!O_x^- + y !^- + z H^+ \rightarrow (\frac{y+1}{2}) !_2 + \frac{z}{2} H_2O$	
			M9
		y =	M10
		y =	,,,,,
	(iv)	Hence, determine the value of x in IO _x ⁻ ion.	
		·	
			!
		x =	M11
, n	45	The maximum error in the volume measured using the pipette is ±0.05 cm ³ .	
(d)	(i)	Calculate the maximum percentage error in the volume of FA 1 used.	
		Calculate the maximum percentage error in the volume of 1771 and 1871	
		Maximum percentage error:	M12

A student suggested that a more accurate value of x could be obtained if a burette is used to measure FA 3 instead of a measuring cylinder. Do you agree with the student? Explain your answer.	
[Total: 13]	M13
	used to measure FA 3 instead of a measuring cylinder. Do you agree with the student? Explain your answer.

2 Determination of the enthalpy change of the decomposition of NaHCO₃

You are required to determine the enthalpy change of the decomposition of NaHCO₃.

$$NaHCO_3 \rightarrow \frac{1}{2}Na_2CO_3 + \frac{1}{2}H_2O + \frac{1}{2}CO_2$$

FA 5 is anhydrous sodium hydrogencarbonate.

FA 6 is anhydrous sodium carbonate.

FA 7 is 1.0 mol dm⁻³ hydrochloric acid.

You will determine the enthalpy change of reaction for each of **FA 5** and **FA 6** with the excess addition of **FA 7**. You will then use these values to calculate the enthalpy change of the decomposition of sodium hydrogencarbonate.

(a) Determining the enthalpy change of reaction between FA 5 and FA 7

Method

- · Weigh an empty weighing boat.
- Weigh about 4.0 g of FA 5 in a weighing boat.
- Using a 50 cm³ measuring cylinder, transfer 50.0 cm³ of FA 7 into a 250 cm³ beaker.
- Stir the solution and record the initial temperature.
- Start timing and do not stop the stopwatch until the whole experiment has been carried out for 7 minutes.
- Record the temperature of the solution mixture every minute for 2 minutes.
- At exactly 3 minutes, carefully add FA 5 into the beaker and stir the mixture. You do not need to measure the temperature of the solution at the 3rd minute. Break up any clumps while stirring.
- Record the temperature at t = 3.5 min.
- Repeat the measurement every 0.5 min until t = 7min.
- Reweigh the weighing boat containing residual sodium hydrogencarbonate.
- In an appropriate format in the space provided below, prepare tables in which to record results for your experiment in (a):
 - o all weighings to an appropriate level of precision,
 - o all values of temperature, T, to an appropriate level of precision,
 - all values of time, t, recorded to the nearest 0.5 min.

Results

(b) Plot a graph of temperature, T, on the y-axis, against time, t, on the x-axis on the grid in Fig. 2.1.

Draw a best-fit straight line taking into account all of the points before t = 3.0 min.

Draw another best-fit straight line taking into account all of the points after the temperature of the mixture has started to rise steadily.

Extrapolate both lines to t = 3.0 min and determine the theoretical fall in temperature.




Fig. 2.1

M16

M17

M18

(c)	Calcu	ulations	
	(i)	Using the answer from (b) , calculate the heat energy taken in during the reaction of FA 5 with FA 7 .	
		[4.2 J is required to raise the temperature of 1 cm³ of solution by 1 °C]	
		-	
		Heat energy taken in =	M20
	(ii)	Calculate the amount of FA 5 , NaHCO ₃ , used in the experiment.	
		[A _r : C, 12.0; H, 1.0; O, 16.0; Na, 23.0]	
		Amount of FA 5 , NaHCO ₃ used =	M21
	(iii)	Hence, determine the enthalpy change of reaction, ΔH_1 for the following reaction.	
	(,	NaHCO ₃ + HC $l \rightarrow$ NaC l + H ₂ O + CO ₂ Δ H ₁	
			M22
		$\Delta H_1 = \dots$	14122

(d) Determine the enthalpy change of reaction between FA 6 and FA 7

Method

- · Weigh an empty weighing boat.
- Place about 2.4 g of FA 6 in the weighing boat.
- Place a dry Styrofoam cup inside a 250 cm³ beaker.
- Using a measuring cylinder, place 50 cm³ of FA 7 into the Styrofoam cup.
- Place the lid onto the cup through the thermometer from the top. Stir the liquid in the cup
 with the thermometer and measure its temperature. Record this temperature in your table.
- Tip cautiously the contents of the weighing boat into the acid in the Styrofoam cup, stir gently with the thermometer and record the highest temperature obtained in your table.
- Calculate the change in temperature.
- · Reweigh the empty weighing boat.
- In an appropriate format in the space provided below, record all measurements of mass and temperature. Include the mass of **FA 6** used in your recording.

Results

M23

M24

M25 M26

M27

(e)	(i)	Given that 4.2 J is required to raise the temperature of 1 cm³ of solution by 1 °C, calculate the heat energy given out in the reaction between FA 6 and FA 7 .	
		Heat energy given out =	M28
	(ii)	Calculate the amount of FA 6 , Na ₂ CO ₃ added to the cup.	
	(ii)	[A _r : C, 12.0; O, 16.0; Na, 23.0]	
		Amount of FA 6 , Na ₂ CO ₃ =	M29
(f)	Calc	ulate the enthalpy change ΔH_2 for the following reaction.	
		$Na_2CO_3 + 2HCI \rightarrow 2NaCI + H_2O + CO_2 \Delta H_2$	
			M 200
		$\Delta H_2 = \dots$	M30 M31
		202	l

	Using your answers in (c)(iii) and (f), calculate the enthalpy change of decomposition of NaHCO ₃ , Δ H ₃ .	
	NaHCO ₃ → ½Na ₂ CO ₃ + ½H ₂ O + ½CO ₂ ΔH ₃	
	ΔH ₃ =	M32
(h)	A student repeats the experiment with 2.0 mol dm ⁻³ HC I instead of the given concentration of 1.0 mol dm ⁻³ , with the same total volume. State and explain the effect of this change on the value of ΔH_3 that she has calculated.	·
	[Total: 20]	M33

3 Qualitative Analysis

You are provided with three solutions, FA 8, FA 9 and FA 10.

Perform the test-tube experiments described below and record your observations in the table.

Use fresh samples of each solution in the tests (a)(i), (a)(iii) and (a)(iv).

		observations with FA 8	observations with FA 9	observations with FA 10
(a)(i)	Place about 1 cm ³ of FA 8 in a test-tube.			
	Add sodium hydroxide, dropwise with shaking, until in excess.			
	Repeat using FA 9 and FA 10. Retain the mixtures for use in (a)(ii).			
(ii)	Add hot water to a 500 cm ³ beaker provided until it is approximately half-filled.			
	Place the test-tubes from (a)(i) in the beaker for 5 minutes.			
	While you are waiting, begin the tests below.	1		
(iii)	Place about 1 cm ³ of FA 8 in a test-tube.			
	Add about 1 cm³ of Fehling's solution and shake the mixture.			
	Place the test-tube in the beaker with hot water for around 5 minutes.			
	While you are waiting, begin test (a)(iv).			
(iv)	Place about 1 cm ³ of FA 8 in a test-tube.			
	Add two drops of acidified potassium manganate(VII) and shake the mixture.			
	Place the test-tube in the beaker with hot water for around 5 minutes.			
	Repeat using FA 9 and FA 10.			

)	 From your observations, state and quote evidence one conclusion you can draw about FA 8 	e for:
	conclusion about FA 8	
	evidence	
	the identity of the cation in FA 9	
	cation in FA 9	
	evidence	
)	Perform the experiment described below and reco	ord your observations in the table.
		observations
	Place about 2 cm³ of FA 10 in a test-tube. Add 2 pieces of magnesium ribbon. Place the test-tube in the beaker with hot water for around 5 minutes.	
t .	Deduce the cation present in FA 10. Hence, so reactions that occurred in (c).	uggest a balanced equation for one of the
	FA 10 contains an anion which is either an iodided Devise and perform a simple test, based on the Opresent in FA 10. Record your observations and 10.	Qualitative Analysis Notes to identify the anion
		[Total: 10]

4 Planning

Lead can form oxides of varied formulae, including Pb_2O and PbO. Lead(II) carbonate, $PbCO_3$, decomposes on heating to form either one of these oxides.

The two possible decomposition equations are:

- 1. $PbCO_3(s) \rightarrow PbO(s) + CO_2(g)$
- 2. $2PbCO_3(s) \rightarrow Pb_2O(s) + 2CO_2(g) + \frac{1}{2}O_2(g)$

You are to plan an experiment to investigate the decomposition of lead(II) carbonate on heating and hence determine the actual decomposition product using gas collection method.

(a) You may assume that you are provided with a sample of lead(II) carbonate and the apparatus commonly found in a college laboratory.

Your plan should include:

- a sketch of the experimental set-up
- calculation of the mass of lead(II) carbonate to be used
- details of how you would carry out the experiment
- a table to show the data you would measure and record during the experiment. Include in your table any other data you would calculate from the experimental results. Insert in your table the letters A, B, C, etc. to represent each data.
- brief, but specific, details of how the results would then be used to reach a conclusion. [A_r: C, 12.0; O, 16.0; Pb, 207.2]

······································

······································	
	M46
	M47
	M48 M49
	M50
	M51

(b)	State one assumption made in the calculations to determine the actual decomposition product.	
		M52
(c)	Suggest another method to collect the gas evolved from the reaction and state its limitation.	
	Method:	
	Limitation:	M53 M54
(d)	Other than the gas collection method, suggest another method to determine the decomposition product.	
		M55
	[Total: 12]	

END OF PAPER

BLANK PAGE

Qualitative Analysis Notes

[ppt. = precipitate]

(a) Reactions of aqueous cations

	react	ion with
cation	NaOH(aq)	NH₃(aq)
aluminium, A/ ³⁺ (aq)	white ppt. soluble in excess	white ppt. insoluble in excess
ammonium, NH₄⁺(aq)	ammonia produced on heating	_
barium, Ba ²⁺ (aq)	no ppt. (if reagents are pure)	no ppt.
calcium, Ca ²⁺ (aq)	white ppt. with high [Ca ²⁺ (aq)]	no ppt.
chromium(III), Cr³+(aq)	grey-green ppt. soluble in excess giving dark green solution	grey-green ppt. insoluble in excess
copper(II), Cu ²⁺ (aq),	pale blue ppt. insoluble in excess	blue ppt. soluble in excess giving dark blue solution
iron(II), Fe ²⁺ (aq)	green ppt. insoluble in excess	green ppt. insoluble in excess
iron(III), Fe ³⁺ (aq)	red-brown ppt. insoluble in excess	red-brown ppt. insoluble in excess
magnesium, Mg²+(aq)	white ppt. insoluble in excess	white ppt. insoluble in excess
manganese(II), Mn ²⁺ (aq)	off-white ppt. insoluble in excess	off-white ppt. insoluble in excess
zinc, Zn²+(aq)	white ppt. soluble in excess	white ppt. soluble in excess

(b) Reactions of anions

ions	reaction	
carbonate, CO ₃ ²⁻	CO₂ liberated by dilute acids	
chloride, CΓ(aq)	gives white ppt. with Ag⁺(aq) (soluble in NH₃(aq))	
bromide, Br ⁻ (aq)	gives pale cream ppt. with Ag ⁺ (aq) (partially soluble in NH ₃ (aq))	
iodide, I ⁻ (aq)	gives yellow ppt. with Ag⁺(aq) (insoluble in NH₃(aq))	
nitrate, NO₃⁻(aq)	NH ₃ liberated on heating with OH ⁻ (aq) and A/ foil	
nitrite, NO₂¯(aq)	NH₃ liberated on heating with OH⁻(aq) and A/ foil; NO liberated by dilute acids (colourless NO → (pale) brown NO₂ in air)	
sulfate, SO ₄ ²⁻ (aq)	gives white ppt. with Ba ²⁺ (aq) (insoluble in excess dilute strong acids)	
sulfite, SO ₃ ²⁻ (aq)	SO ₂ liberated with dilute acids; gives white ppt. with Ba ²⁺ (aq) (soluble in dilute strong acids)	

(c) Test for gases

ions	reaction		
ammonia, NH₃	turns damp red litmus paper blue		
carbon dioxide, CO ₂	gives a white ppt. with limewater (ppt. dissolves with excess CO ₂)		
chlorine, C/2	bleaches damp litmus paper		
hydrogen, H ₂	"pops" with a lighted splint		
oxygen, O ₂	relights a glowing splint		
sulfur dioxide, SO ₂	turns aqueous acidified potassium manganate(VII) from purple to colourless		

(d) Colour of halogens

halogen	colour of element	colour in aqueous solution	colour in hexane	
chlorine, Cl ₂	greenish yellow gas	pale yellow	pale yellow	
bromine, Br ₂	reddish brown gas / liquid	orange	orange-red	
iodine, I ₂	black solid / purple gas	brown	purple	

			Clas	ss Adm No
Candidate Name: _	_			
·				Shift
	-11	•		

N	millennia institute		
---	------------------------	--	--

Shift	
Laboratory	
	ŀ

2021 Preliminary ExaminationsPre-University 3

H2 CHEMISTRY

9729/04

Paper 4 Practical

1 Sept 2021

2 hours 30 min

Candidates answer on the Question paper.

READ THESE INSTRUCTIONS FIRST

Do not turn over this question paper until you are told to do so

Write your name, class and admission number on all the work you hand in.

Write in dark blue or black pen.

You may use an HB pencil for any diagrams or graphs.

Do not use staples, paper clips, glue or correction fluid.

Answer all questions in the spaces provided on the Question Paper.

The use of an approved scientific calculator is expected, where appropriate.

You may lose marks if you do not show your working or if you do not use appropriate units.

Qualitative Analysis Notes are printed at the back of the Question Paper.

At the end of the examination, fasten all your work securely together.

The number of marks is given in brackets [] at the end of each question or part question.

Question	1	2	3	4	Total
Marks					
	13	20	10	12	55

1 Determine the formula of the ion IO_x by titration

In this experiment you will determine the formula of the ion, IOx-.

lodine, l₂, is formed when IO_x ions reacts with an excess of iodide ions.

The equation for this reaction is:

$$IO_x^- + y I^- + z H^+ \rightarrow (\frac{y+1}{2}) I_2 + \frac{z}{2} H_2O$$

where x, y and z are all integers.

The amount of iodine produced will then be determined by titration with thiosulfate ions, S₂O₃²⁻.

$$l_2 + 2S_2O_3^{2-} \rightarrow 2I^- + S_4O_6^{2-}$$

FA 1 is a solution containing 0.0150 mol $dm^{-3} IO_x^-$ ions.

FA 2 is dilute sulfuric acid, H₂SO₄.

FA 3 is 0.500 mol dm⁻³ potassium iodide, KI.

FA 4 is $0.100 \text{ mol dm}^{-3}$ sodium thiosulfate, $Na_2S_2O_3$.

starch indicator

(a) Method

- Pipette 25.0 cm³ of FA 1 into a conical flask.
- Using a measuring cylinder, add 20.0 cm³ of FA 2 to the conical flask.
- Using a measuring cylinder, add 10.0 cm³ of FA 3 to the conical flask. The solution will turn brown as iodine is produced.
- Fill the burette with FA 4.
- Add FA 4 from the burette until the solution in the conical flask turns yellow.
- Add 1 cm³ of starch indicator to the conical flask. The solution will turn blue-black.
- Continue to add more FA 4 dropwise from the burette until the blue-black colour just disappears.

This is the end-point of the titration.

- Carry out as many accurate titrations as you think necessary to obtain consistent results.
- Record the burette readings and the volume of FA 4 added in each accurate titration in the space below.

Results

Initial burette reading / cm ³	0.00	0.00
Final burette reading / cm ³	22.50	22.50
Volume of FA 4 added / cm ³	22.50	22.50
70.000	√	✓

M1 table w headers and units

m2 readings to 0.05 cm³ + correct calc

M3
2
consistent
values
within
0.10

M4 accuracy

M5 accuracy (b) From your titration results, obtain a value for the volume of FA 4 to be used in your calculations. Show clearly how you obtained this value.

Volume of **FA 4** added =
$$\frac{22.50+22.50}{2}$$

= 22.50 cm³

M6 correct calc w working

25.0 cm³ of **FA 1** required <u>22.50</u> cm³ of **FA 4**.

- (c) Calculations
 - (i) Use your answer to (b) to calculate the amount of iodine that was formed.

Amount of
$$S_2O_3^{2-} = 0.100 \times \frac{22.50}{1000}$$

= 0.00225 mol
Amount of $I_2 = \frac{1}{2} \times 0.00225$
= 1.13 x 10⁻³ mol (3sf)

Amount of $I_2 = 1.13 \times 10^{-3} \text{ mol}$

М7

(ii) Calculate the amount of IO_x^- ions in 25.0 cm³ of FA 1.

Amount of
$$IO_x^- = 0.0150 \times \frac{25.0}{1000}$$

= 3.75 x 10⁻⁴ mol

amount of IO_x^- ions= 3.75×10^{-4} mol

M8

(iii) Using your answers in (c)(i) and (c)(ii), calculate the value of y.

$$\text{IO}_x^- + \textbf{\textit{y}} \text{ I}^- + \textbf{\textit{z}} \text{ H}^+ \rightarrow \big(\frac{\textbf{\textit{y}+1}}{2}\big) \text{ I}_2 + \frac{\textbf{\textit{z}}}{2} \text{ H}_2 \text{O}$$

$$\frac{\text{Amount of I}_2}{\text{Amount of IO}_{\overline{x}}} = \frac{y+1}{2}$$

$$\frac{1.125 \times 10^{-3}}{3.75 \times 10^{-4}} = \frac{y+1}{2}$$

y = 5 (must round off to nearest integer)

М9

M10 (c) ans to 3sf / integer + correct

units

y = 5

(iv) Hence, determine the value of x in IO_x^- ion.

$$10x^{-} + 51^{-} + z H^{+} \rightarrow 31_{2} + \frac{z}{2} H_{2}O$$

Total charge on LHS = Total charge on RHS

$$-1 - 5 + z = 0 + 0$$

z = 6

$$10x^{-} + 51^{-} + 6H^{+} \rightarrow 31z + 3H_{2}O$$

x = 3

x = <u>3</u>

M11

(d) (i) The maximum error in the volume measured using the pipette is ±0.05 cm³.Calculate the maximum percentage error in the volume of FA 1 used.

Max % error =
$$\frac{0.05}{25.0}$$
 x 100% = 0.200 %

Maximum percentage error: $\underline{0.200~\%}$

M12

(ii) A student suggested that a more accurate value of **x** could be obtained if a burette is used to measure **FA 3** instead of a measuring cylinder.

Do you agree with the student? Explain your answer.

Disagree. FA 3 is added in excess, hence there is no need to use a precise instrument.

M13

[Total: 13]

2 Determination of the enthalpy change of the decomposition of NaHCO₃

You are required to determine the enthalpy change of the decomposition of NaHCO₃.

$$NaHCO_3 \rightarrow \frac{1}{2}Na_2CO_3 + \frac{1}{2}H_2O + \frac{1}{2}CO_2$$

FA 5 is anhydrous sodium hydrogencarbonate.

FA 6 is anhydrous sodium carbonate.

FA 7 is 1.0 mol dm⁻³ hydrochloric acid.

You will determine the enthalpy change of reaction for each of **FA 5** and **FA 6** with the excess addition of **FA 7**. You will then use these values to calculate the enthalpy change of the decomposition of sodium hydrogencarbonate.

(a) Determining the enthalpy change of reaction between FA 5 and FA 7

 $NaHCO_3 + HC/ \rightarrow NaC/ + H_2O + CO_2$

Method

- · Weigh an empty weighing boat.
- Weigh about 4.0 g of FA 5 in a weighing boat.
- Using a 50 cm³ measuring cylinder, transfer 50.0 cm³ of FA 7 into a 250 cm³ beaker.
- Stir the solution and record the initial temperature.
- Start timing and do not stop the stopwatch until the whole experiment has been carried out for 7 minutes.
- Record the temperature of the solution mixture every minute for 2 minutes.
- At exactly 3 minutes, carefully add FA 5 into the beaker and stir the mixture. You do not need to measure the temperature of the solution at the 3rd minute. Break up any clumps while stirring.
- Record the temperature at t = 3.5 min.
- Repeat the measurement every 0.5 min until t = 7min.
- Reweigh the weighing boat containing residual sodium hydrogencarbonate.
- In an appropriate format in the space provided below, prepare tables in which to record results for your experiment in (a):
 - o all weighings to an appropriate level of precision,
 - o all values of temperature, T, to an appropriate level of precision,
 - o all values of time, t, recorded to the nearest 0.5 min.

Results

Mass of empty weighing boat / g	3.11
Mass of weighing boat with FA 5 / g	7.11
Mass of weighing boat with residual FA 5 / g	3.13
Mass of FA 5 transferred / g	3.98
Mass of A 3 fransierred 9	

į	t / min	0.0	10	20	3.5	4.0	4.5	5.0	5.5	6.0	6.5	7.0
1	T/°C	0.0	1.0	2.0	0.0	25.0	05.4	25.4	25.5	25.6	25.8	25.8
1	T/°C	30.0	30.0	30.0	27.2	25.6_	25.4	25.4	25.5	23.0	20.0	20.0

M14 3 mass readings + lime, temp


> M15 units

(b) Plot a graph of temperature, T, on the y-axis, against time, t, on the x-axis on the grid in Fig. 2.1.

Draw a best-fit straight line taking into account all of the points before t = 3.0 min.

Draw another best-fit straight line taking into account all of the points after the temperature of the mixture has started to rise steadily.

Extrapolate both lines to t = 3.0 min and determine the theoretical fall in temperature.



M16 axes labelled

> M17 scale

with units

M18 all points plotted correctly

M19
2 best fit lines, extrapolated to t = 3.0, 2nd line must only use increasing / same

value ots

- (c) Calculations
 - (i) Using the answer from (b), calculate the heat energy taken in during the reaction of FA 5 with FA 7.

[4.2 J is required to raise the temperature of 1 cm³ of solution by 1 °C]

Q =
$$mc\Delta T$$

= $(50.0)(4.2)(5.0)$
= $1050 J$

Heat energy taken in = 1050 J

M20

(ii) Calculate the amount of FA 5, NaHCO₃, used in the experiment.

[A_r: C, 12.0; H, 1.0; O, 16.0; Na, 23.0]

Amount of NaHCO₃ =
$$\frac{3.98}{23.0+1.0+12.0+3(16.0)}$$

= 0.0474 mol

Amount of FA 5, NaHCO₃ used = 0.0474 mol

M21

(iii) Hence, determine the enthalpy change of reaction, ΔH_1 for the following reaction. NaHCO₃ + HCl \rightarrow NaCl + H₂O + CO₂ ---- Δ H₁

$$\Delta H_1 = \frac{Q}{n_{LR}} \times CLR$$

= $\frac{1050}{0.04738} \times 1$
= $+22.2 \text{ kJ mol}^{-1}$

 $\Delta H_1 = +22.2 \text{ kJ mol}^{-1}$

M22

(d) Determine the enthalpy change of reaction between FA 6 and FA 7

Method

- · Weigh an empty weighing boat.
- Place about 2.4 g of FA 6 in the weighing boat.
- Place a dry Styrofoam cup inside a 250 cm³ beaker.
- Using a measuring cylinder, place 50 cm³ of FA 7 into the Styrofoam cup.
- Place the lid onto the cup through the thermometer from the top. Stir the liquid in the cup
 with the thermometer and measure its temperature. Record this temperature in your
 table.
- Tip cautiously the contents of the weighing boat into the acid in the Styrofoam cup, stir gently with the thermometer and record the highest temperature obtained in your table.
- Calculate the change in temperature.
- · Reweigh the empty weighing boat.
- In an appropriate format in the space provided below, record all measurements of mass and temperature. Include the mass of **FA 6** used in your recording.

Results

Mass of empty weighing boat / g	2.85
Mass of weighing boat with FA 6 / g	5.25
Mass of weighing boat with residual FA 6 / g	2.87
Mass of FA 6 transferred / g	2.38

Initial temperature / °C	30.0
Final temperature / °C	33.0
Change in temperature / °C	3.0

M23

3 mass
readings +
1 mass
calc + 2
temp
readings +
1 temp
calc

M24

2.35 ≤
mass
transferr
ed ≤
2.45 +
correct
calcs

M25

all mass readings 2dp, all temp readings to 0.1 °C

M26 accuracy

M27 accuracy

(e) (i) Given that 4.2 J is required to raise the temperature of 1 cm³ of solution by 1 °C, calculate the heat energy given out in the reaction between FA 6 and FA 7.

$$Q = mc\Delta T$$
= (50.0)(4.2)(3.0)
= 630 J

Heat energy given out = 630 J

M28

(ii) Calculate the amount of FA 6, Na₂CO₃ added to the cup.

[Ar: C, 12.0; O, 16.0; Na, 23.0]

Amount of NaHCO₃ =
$$\frac{2.38}{2(23.0)+12.0+3(16.0)}$$

= 0.0225 mol

Amount of FA 6, Na₂CO₃ = 0.0225 mol

M29

(f) Calculate the enthalpy change ΔH_2 for the following reaction.

$$Na_2CO_3 + 2HCl \rightarrow 2NaCl + H_2O + CO_2 -\!\!\!\!-\!\!\!\!- \Delta H_2$$

$$\Delta H_1 = \frac{Q}{n_{LR}} \times CLR$$
$$= \frac{630}{0.0225} \times 1$$
$$= -28.1 \text{ kJ mol}^{-1}$$

M30

M31 correct signs for (c)(iii) and (f)

 $\Delta H_2 = -28.1 \text{ kJ mol}^{-1}$

(g) Using your answers in (c)(iii) and (f), calculate the enthalpy change of decomposition of NaHCO₃, ΔH₃.

$$NaHCO_3 \rightarrow \frac{1}{2}Na_2CO_3 + \frac{1}{2}H_2O + \frac{1}{2}CO_2 ---- \Delta H_3$$

HC/ + NaHCO₃
$$\rightarrow$$
 ½Na₂CO₃ + ½H₂O + ½CO₂ + HC/
$$\Delta H_1 \qquad \qquad \Delta H_2$$
 NaC/ + H₂O + CO₂

$$\Delta H_3 = \Delta H_1 - \frac{1}{2}\Delta H_2$$

= +22.2 - $\frac{1}{2}$ (-28.1)
= +36.2 kJ mol⁻¹

 $\Delta H_3 = +36.2 \text{ kJ mol}^{-1}$

M32

(h) A student repeats the experiment with 2.0 mol dm⁻³ HC*l* instead of the given concentration of 1.0 mol dm⁻³, with the same total volume.

State and explain the effect of this change on the value of ΔH_3 that she has calculated.

No effect. This is because \underline{HCI} was already used in excess, hence $\underline{\Delta H_1}$ and $\underline{\Delta H_2}$ calculated is independent of the amount of \underline{HCI} . Therefore, $\underline{\Delta H_3}$ will not be affected.

M33

[Total: 20]

3 Qualitative Analysis

You are provided with three solutions, FA 8, FA 9 and FA 10.

Perform the test-tube experiments described below and record your observations in the table.

Use fresh samples of each solution in the tests (a)(i), (a)(iii) and (a)(iv).

		observations with FA 8	observations with FA 9	observations with FA 10
(a)(i)	Place about 1 cm ³ of FA 8 in a test-tube. Add sodium hydroxide, dropwise with shaking, until in excess.	no ppt (1)	off-white / pale brown / beige / buff ppt and insoluble in	green ppt and insoluble in excess / ppt darkens (1)
:	Repeat using FA 9 and FA 10. Retain the mixtures for use in (a)(ii).		excess / ppt darkens (1)	
(ii)	Add hot water to a 500 cm ³ beaker provided until it is approximately half-filled.			brown / black ppt (1) [ppt must be darker
	Place the test-tubes from (a)(i) in the beaker for 5 minutes.	yellow orange / brown solution (1)	brown ppt (1)	than previously] gas turns red litmus
	While you are waiting, begin the tests below.			blue (1)
(iii)	Place about 1 cm ³ of FA 8 in a test-tube.			
	Add about 1 cm ³ of Fehling's solution and shake the mixture.	brick-red / red /		
	Place the test-tube in the beaker with hot water for around 5 minutes.	orange ppt (1)		
	While you are waiting, begin test (a)(iv).			
(iv)	Place about 1 cm ³ of FA 8 in a test-tube.			
	Add two drops of acidified potassium manganate(VII) and shake the mixture.	purple to colourless solution / brown ppt	brown ppt (1)	purple to colourless solution (1)
	Place the test-tube in the beaker with hot water for around 5 minutes.	(1)		
	Repeat using FA 9 and FA 10.			

- (b) From your observations, state and quote evidence for:
 - one conclusion you can draw about FA 8

conclusion about FA 8 contains aliphatic aldehyde / is a reducing agent

evidence brick-red ppt forms when Fehling's solution was added

the identity of the cation in FA 9

cation in FA 9 Mn2+

evidence off-white ppt upon addition of NaOH, which is insoluble in excess NaOH

M38

M39

(c) Perform the experiment described below and record your observations in the table.

	observations	
 Place about 2 cm ³ of FA 10 in a test-tube.	Gas evolved causes lighted splint to 'pop'	
Add 2 pieces of magnesium ribbon.	OR	
Place the test-tube in the beaker with hot water		N
for around 5 minutes.	Green ppt formed / Brown/black deposit on Mg / brown or black ppt	14

M40

(d) Deduce the cation present in **FA 10**. Hence, suggest a balanced equation for one of the reactions that occurred in (c).

FA 10 contains Fe2+

M41

$$2Fe(OH)_2 + [O] + H_2O \rightarrow 2Fe(OH)_3$$
 OR

$$Fe^{2+} + Mg \rightarrow Fe + Mg^{2+}$$

(e) FA 10 contains an anion which is either an iodide ion or sulfate ion.

Devise and perform a simple test, based on the Qualitative Analysis Notes to identify the anion present in **FA 10**. Record your observations and hence deduce the identity of the anion in **FA 10**.

To a test-tube containing 1 cm 3 **FA 10**, add 1 cm 3 of BaC l_2 solution. A white precipitate is formed, hence the anion is sulfate ion.

OR

To a test-tube containing 1 cm³ **FA 10**, add 1 cm³ of AgNO₃ solution, followed by excess NH₃. No yellow precipitate is formed, hence the anion is sulfate ion.

M42 M43

[Total: 10]

4 Planning

Lead can form oxides of varied formulae, including Pb_2O and PbO. Lead(II) carbonate, $PbCO_3$, decomposes on heating to form either one of these oxides.

The two possible decomposition equations are:

1.
$$PbCO_3(s) \rightarrow PbO(s) + CO_2(g)$$

2.
$$2PbCO_3(s) \rightarrow Pb_2O(s) + 2CO_2(g) + \frac{1}{2}O_2(g)$$


You are to plan an experiment to investigate the decomposition of lead(II) carbonate on heating and hence determine the actual decomposition product using gas collection method.

(a) You may assume that you are provided with a sample of lead(II) carbonate and the apparatus commonly found in a college laboratory.

Your plan should include:

- a sketch of the experimental set-up
- calculation of the mass of lead(II) carbonate to be used
- details of how you would carry out the experiment
- a table to show the data you would measure and record during the experiment. Include
 in your table any other data you would calculate from the experimental results. Insert
 in your table the letters A, B, C, etc. to represent each data.
- brief, but specific, details of how the results would then be used to reach a conclusion. $[A_r:C,\,12.0;\,O,\,16.0;\,Pb,\,207.2]$

Apparatus:



M44

Apparatus:
Sketch of
setup
generally
ok (either
gas
syringe or
displacem
ent of
water)

M45

Apparatus: Key apparatus labelled (Bunsen burner. volume. rictionless gas syringe OR measuring cylinder / burette with water trough

Pre-calculations:

Since 100 cm³ gas syringe is used, assume that 80 cm³ of gas is collected.

Amount of gas =
$$\frac{80}{24000}$$
 = 0.00333 mol

If decomposition via equation 1,

Amount of PbCO₃ = Amount of gas (CO₂ only)

$$= 0.00333 \text{ mol}$$

Mass of PbCO₃ required = $0.00333 \times [207.2+12.0+3(16.0)]$

$$\bar{}$$
 = 0.891 g

If decomposition via equation 2,

Amount of PbCO₃ = $\frac{2}{2.5}$ x Amount of gas (CO₂ and O₂)

$$= 0.00267 \text{ mol}$$

Mass of PbCO₃ required = $0.00267 \times [207.2+12.0+3(16.0)]$

$$= 0.713 q$$

To ensure that volume does not exceed 80 cm³, use smaller mass of 0.713 g.

Procedure:

- 1. Using a weighing balance, weigh 0.71 g of PbCO₃ and transfer it into a boiling tube.
- 2. Set up apparatus as shown in diagram on previous page.
- 3. Record initial volume of the gas syringe in a table.
- 4. <u>Light the Bunsen burner to begin heating. Turn off the Bunsen burner when the plunger of the gas syringe no longer moves.</u>
- 5. Let the set-up cool, then record final volume of the gas syringe.

Table of measurements:

Initial volume of gas syringe / cm ³	Α
Final volume of gas syringe / cm ³	В
Volume of gas evolved / cm ³	C = B - A

Conclusions:

Amount of PbCO₃ =
$$\frac{0.71}{207.2+12.0+3(16.0)}$$
 = 0.00266 mol

Amount of gas =
$$\frac{c}{24000}$$
 mol

Ratio of amount of gas to PbCO₃ = $\frac{C}{63.77}$

If ratio is closer to 1, product is PbO. If ratio is closer to 1.25, product is PbO₂.

M46

Pre-calc: Ensures choice of mass is linked to volume of apparatus (70-80%)

M47

Pre-calc:
Determine
s that
smaller
mass of
solid (from
eqn 2) is
more
appropriat
e

M48

Procedure: Weighs and heats

M49

Procedure:
Recording
of initial
volume
and final
volume
only when
no further
change

M50

Recording: Table of volume measurem ents + calc with units and data letters

M51

Conclusion
: Calc for
ratio of
amount of
gas to
amount of
PbCO₃
used

(b) State one assumption made in the calculations to determine the actual decomposition product.

Gas syringe: The volume of gas collected is at room temperature and pressure Water displacement: No gas dissolved in the water

M52

(c) Suggest another method to collect the gas evolved from the reaction and state its limitation.

Water displacement: CO₂ is slightly soluble in water (1), hence some gas may have been lost (1) which makes the data collected unreliable.

M53

OR

M54

Gas syringe: There is some friction in the gas syringe (1), which results in the plunger not moving properly, (1) hence making the data collected unreliable.

(d) Other than the gas collection method, suggest another method to determine the decomposition product.

Gravimetry / test the gas with a glowing splint.

M55

[Total: 12]

END OF PAPER

BLANK PAGE

Qualitative Analysis Notes

[ppt. = precipitate]

(a) Reactions of aqueous cations

	reaction with				
cation	NaOH(aq)	NH₃(aq)			
aluminium, A/³⁺(aq)	white ppt. soluble in excess	white ppt. insoluble in excess			
ammonium, NH₄⁺(aq)	ammonia produced on heating	-			
barium, Ba ²⁺ (aq)	no ppt. (if reagents are pure)	no ppt.			
calcium, Ca ²⁺ (aq)	white ppt. with high [Ca ²⁺ (aq)]	no ppt.			
chromium(III), Cr ³⁺ (aq)	grey-green ppt. soluble in excess giving dark green solution	grey-green ppt. insoluble in excess			
copper(II), Cu ²⁺ (aq),	pale blue ppt. insoluble in excess	blue ppt. soluble in excess giving dark blue solution			
iron(II), Fe ²⁺ (aq)	green ppt. insoluble in excess	green ppt. insoluble in excess			
iron(III), Fe ³⁺ (aq)	red-brown ppt. insoluble in excess	red-brown ppt. insoluble in excess			
magnesium, Mg ²⁺ (aq)	white ppt. insoluble in excess	white ppt. insoluble in excess			
manganese(II), Mn ²⁺ (aq)	off-white ppt. insoluble in excess	off-white ppt. insoluble in excess			
zinc, Zn ²⁺ (aq)	white ppt. soluble in excess	white ppt. soluble in excess			

(b) Reactions of anions

ions	reaction		
carbonate, CO ₃ ²⁻	CO₂ liberated by dilute acids		
chloride, C/⁻(aq)	gives white ppt. with Ag⁺(aq) (soluble in NH₃(aq))		
bromide, Br ⁻ (aq)	gives pale cream ppt. with Ag⁺(aq) (partially soluble in NH₃(aq))		
iodide, ⊺(aq)	gives yellow ppt. with Ag⁺(aq) (insoluble in NH₃(aq))		
nitrate, NO ₃ -(aq)	NH ₃ liberated on heating with OH ⁻ (aq) and A/ foil		
nitrite, NO ₂ (aq)	NH₃ liberated on heating with OH⁻(aq) and A/ foil; NO liberated by dilute acids (colourless NO → (pale) brown NO₂ in air)		
sulfate, SO ₄ ²⁻ (aq)	gives white ppt. with Ba ²⁺ (aq) (insoluble in excess dilute strong acids)		
sulfite, SO ₃ ²⁻ (aq)	SO ₂ liberated with dilute acids; gives white ppt. with Ba ²⁺ (aq) (soluble in dilute strong acids)		

(c) Test for gases

ions	reaction			
ammonia, NH ₃	turns damp red litmus paper blue			
carbon dioxide, CO ₂	gives a white ppt. with limewater (ppt. dissolves with excess CO ₂)			
chlorine, CI ₂	bleaches damp litmus paper			
hydrogen, H₂	"pops" with a lighted splint			
oxygen, O ₂	relights a glowing splint			
sulfur dioxide, SO ₂	turns aqueous acidified potassium manganate(VII) from purple to colourless			

(d) Colour of halogens

halogen	colour of element	colour in aqueous solution	colour in hexane	
chlorine, Cl ₂	greenish yellow gas	pale yellow	pale yellow	
bromine, Br ₂	reddish brown gas / liquid	orange	orange-red	
iodine, I2	black solid / purple gas	brown	purple	

www.testpapersfree.com